bio nova

Omniflow™ II

Instructions for Use / Gebrauchsanleitung / Gebruiksaanwijzing /
Mode d'emploi / Istruzioni d'uso / Instrucciones de uso / Instruções de utilização
English – German – Dutch – French – Italian – Spanish – Portuguese

CE 0197

BN 57-3(5)-3

Omniflow™ II

Contents			PAGE
Vascular prosthesis	Instructions for Use	English	2
Gefäßprothese	Gebrauchsanweisung	Deutsch	8
Vaatprothese	Gebruiksaanwijzing	Nederlands	14
Prothèse vasculaire	Mode d'Emploi	Français	20
Protesi vascolare	Istruzioni d'uso	Italiano	26
Prótesis vascular	Instrucciones de uso	Español	32
Prótese vascular	Instruções de utilização	Português	38

bio nova

STERILE – FOR SINGLE USE ONLY

Description

Omniflow II Vascular Prosthesis is constructed from a composite biosynthetic material comprised of polyester mesh incorporated within a cross-linked ovine fibrocollagenous tissue matrix. The prosthesis is sterilised in a glutaraldehyde solution.

The prosthesis is supplied sterile and non-pyrogenic in a solution of 50% ethanol. The prosthesis remains sterile unless the primary package is opened or damaged.

The **Omniflow II Straight Vascular Prosthesis** is mounted on a glass mandrel contained in a glass tube. The mandrel design prevents the prosthesis slipping off the mandrel when it is removed from the glass tube. The diameter and minimum length of the prosthesis is specified on the label applied to the glass tube.

The **Omniflow II Curved Vascular Prosthesis** is contained in a sterile flexible inner bag within an outer bag. The diameter and minimum length of the prosthesis is specified on the label applied to the outer surface of the outer bag.

Indications for Use

The Omniflow II vascular prosthesis is indicated for the replacement, reconstruction, bypassing or patching of diseased vessels in patients suffering occlusive or aneurismal diseases, in trauma patients requiring vascular replacement or for patients requiring vascular access such as for haemodialysis.

The **Omniflow II Straight Vascular Prosthesis** is intended:

1. To bypass, replace or reconstruct diseased or injured blood vessels
2. To patch and repair peripheral vessels
3. For arteriovenous access when a straight configuration is required.

The **Omniflow II Curved Vascular Prosthesis** is intended for arteriovenous access when a looped configuration is required.

Storage

Storage recommended at room temperature.

Warnings

1. Do NOT re-sterilise the Omniflow II prosthesis. It is supplied sterile and pyrogen free. Use the prosthesis immediately after opening the package and discard any unused portions.
2. Do NOT expose the prosthesis to sunlight or extreme temperatures (Do not freeze, do not store in temperatures above 55°C).
3. Do NOT use the prosthesis if the primary package is damaged as sterility may be compromised.
4. Do NOT use the prosthesis if the glass mandrel is broken.
5. Do NOT use the prosthesis if it is not completely covered by the storage solution.
6. Do NOT attempt to reposition the prosthesis after removal of the tunnelling instrument.
7. Do NOT straighten the curved prosthesis during preparation or implantation, as this will cause disruption of the mesh tissue interface.
8. Do NOT use the straight prosthesis to fashion a looped arteriovenous access as this may cause kinking.
9. Do NOT pull, stretch, twist, squeeze or pinch the body of the prosthesis.
10. Do NOT use ablation techniques such as cutting balloons, laser, or radio frequency ablation with the Omniflow II prosthesis.
11. Do NOT attempt to dilate the prosthesis with balloon angioplasty or stenting procedures.

Technical Information / Precautions

1. Ensure the rinsing procedure has been performed to remove the storage solution prior to implanting the prosthesis. Failure to do so may cause occlusion. Keep the prosthesis moist with sterile physiological saline during the procedure.
2. The use of a tunnelling instrument for the passage of the prosthesis is essential. Failure to do so may cause disruption to the biosynthetic material and lead to occlusion, dilatation or aneurysm formation.
3. Ensure that the prosthesis does not become twisted when passing through the tunnelling instrument as this may lead to occlusion.
4. Avoid cross clamping with metal instruments as this may damage the prosthesis and cause occlusion, dilatation or aneurysm formation. If clamping is necessary use only a-traumatic clamps and avoid repeated or excessive clamping in the same position on the prosthesis.
5. The prosthesis has minimal longitudinal elasticity. Ensure the prosthesis is cut to the correct length. If it is too short it may cause suture pullout with a risk of anastomotic aneurysm. If it is too long it may kink and cause occlusion.
6. Cut off the sections of the prosthesis which were clamped during rinsing. Ensure that the full wall thickness and a mesh eyelet are incorporated with each stitch when performing the anastomosis. Failure to do so may result in stitch pullout and anastomotic aneurysm formation.
7. When the prosthesis is used for arteriovenous access some redness and swelling may be present over the implant area for a few days following implant.

8. Insufficient data is available on which to base any conclusions regarding the use of the Omniflow II vascular prosthesis for aortocoronary bypass procedures.
9. It is recommended that the Omniflow II prosthesis not be implanted into an area with an ongoing infection.

Thrombectomy

Prompt intervention when occlusion occurs has been found to restore prosthetic function. Therefore the patient should be advised to report to the clinician if symptoms return.

When performing a thrombectomy choose an embolectomy catheter of the appropriate size and follow the catheter manufacturer's instructions for use.

To avoid damage to the wall and flow surface of the prosthesis:

1. Do NOT exert undue pressure when inflating the balloon
2. Do NOT exert undue force when withdrawing the catheter from the prosthesis.

Contraindications

The prosthesis should not be used in patients with a known hypersensitivity to ovine material or glutaraldehyde.

Potential Complications with the use of Vascular Prostheses

Complications may occur with the use or in conjunction with any vascular prosthesis and include but are not limited to: infection, thrombosis, dilatation, aneurysm formation, pseudoaneurysm formation, suture pullout and adverse tissue responses.

Late aneurysm formation in Omniflow II has been reported (more than 4 years after implantation). Ongoing patient monitoring is recommended.

Wall integrity of the prosthesis may be adversely affected by collagenase-producing microorganisms.

Pre-implant preparation

Prepare a sterile basin with adequate sterile physiological saline to cover the prosthesis, a sterile 20 ml syringe, sterile heparin 5,000 I.U. per ml.

Removing the prostheses from the container

STRAIGHT VASCULAR PROSTHESIS

1. Remove the seal and cap from the glass tube. Remove the silicone stopper with sterile forceps. The forceps and stopper are discarded. Fig. 1.
2. Aseptically remove the prosthesis from the glass tube by grasping the end of the glass mandrel with sterile forceps and lifting it gently out of the glass tube. Fig. 2. The hook on the bottom of the mandrel will prevent the graft from slipping off the mandrel.
3. Grasp the mandrel by the hook, turn it upside down and allow the prosthesis to slide into the basin of sterile saline. If the prosthesis does not slide off easily, gently push it down the mandrel. Do NOT pull the prosthesis.

CURVED VASCULAR PROSTHESIS

1. Aseptically remove the sterile inner bag containing the prosthesis from the outer bag. Fig. 3. Only the inner bag is sterile.

2. Cut the corner of the inner bag with sterile scissors and allow the storage solution to drain into a container. Fig. 4. Discard the container of storage solution.
3. Remove the prosthesis from the inner bag and place it into the basin of sterile physiological saline.

Rinsing the prosthesis

STRAIGHT VASCULAR PROSTHESIS

1. Use a sterile syringe to flush the lumen of the prosthesis with 20 ml of sterile physiological saline. Fig. 5. Repeat at least twice.
2. Clamp one end of the prosthesis. Insert sterile heparin at a concentration of 5,000 I.U. per ml into the lumen of the prosthesis to a maximum of 50,000 I.U. (10 ml) for prostheses 35 cm and longer and to a maximum of 25,000 I.U. (5 ml) for the shorter prostheses. Add sufficient sterile saline if necessary to fill the prosthesis so that the entire flow surface of the graft is coated with the solution.
3. Apply a clamp to the other end of the prosthesis.
4. Place the prosthesis back into the basin completely covered by the sterile saline until it is time for implantation. This should be a minimum of 10 minutes.

Fig. 1

Fig. 2

Fig. 3


Fig. 4

When it is time to implant the prosthesis, remove the clamps from the prosthesis and allow the heparin and saline to drain out.

CURVED VASCULAR PROSTHESIS

The curved prosthesis is rinsed in the same manner as described above for the straight prosthesis. The curved configuration is maintained throughout the procedure. Fig. 6.

Do NOT straighten the prosthesis.

Implantation

Prepare the incisions in the normal fashion.

Use of a hollow tunnelling instrument

STRAIGHT VASCULAR PROSTHESIS

Form a tunnel for the prosthesis using a metal or plastic hollow tunnelling instrument and pass the prosthesis through the instrument. Fig. 7.

CURVED VASCULAR PROSTHESIS

Marking the proposed graft layout on the skin with a sterile surgical marker may facilitate final positioning of the prosthesis. Maintain the curve of the graft by inserting 2 hollow tunnelling instruments and passing both straight sections of the graft through simultaneously from the distal to the proximal incision. Fig. 8. Ensure that the graft does not kink at the curve during placement.

Carefully remove the tunnelling instrument when the graft has been correctly placed.

Anastomoses

1. When fashioning the ends for the anastomoses it is essential to cut off those portions of the prosthesis that were clamped during preparation.
2. Monofilament polypropylene is the preferred suture material. During suturing evert the edges of the prosthesis and ensure the full wall thickness and a mesh eyelet are taken up with each stitch. Fig. 9.
3. Prior to completion of the second anastomosis flush the prosthesis with blood to ensure complete removal of the heparin.
4. An intra-operative angiogram may be performed at the time of the procedure to document function. Ensure that the artery rather than the prosthesis is used for the injection.

At the completion of the procedure the patient record labels supplied with the prosthesis should be affixed in the patient and hospital records.

Fig. 5

Fig. 6


Fig. 7


Fig. 8


Fig. 9

Arteriovenous Access

1. It is preferable to allow the prosthesis to heal for a minimum of two weeks prior to commencement of puncturing.
2. Do NOT puncture the arteriovenous access repeatedly in the same site as this may lead to disruption of the wall of the prosthesis, haematoma formation or pseudoaneurysm formation. Rotation of the puncture sites along the length of the prosthesis is necessary.
3. Do NOT puncture the curve of the prosthesis.
4. Do NOT puncture within 3 cm of the anastomotic sites.
5. Strict adherence to aseptic technique is essential to minimise the risk of infection.
6. Moderate compression of the puncture site will achieve haemostasis.

Vascular Patch

Cut the tubular prosthesis longitudinally to fashion into a vascular patch. Ensure that the blood flow surface is not damaged. When suturing the patch into place, ensure that the full wall thickness and the mesh are incorporated into each stitch.

Return to Bio Nova International Pty Ltd

Patient Implant Data Form

The implant data form should be completed immediately after surgery by the surgeon and returned in the envelope provided.

Explanted specimens of prosthesis

If the Omniflow II vascular prosthesis is removed from a patient for any reason it may be returned to Bio Nova International Pty Ltd for analysis. Rinse blood and contaminants from the prosthesis using sterile saline. Return the specimen in a sealed container of 10% formalin. Include the serial number, catalogue number and reason for removal. Contact your local Omniflow II representative who will make arrangements for the specimen to be returned.

	Attention, see instructions for use	REF	Catalogue number
STERILE	Chemically sterilised in glutaraldehyde. Contents sterile unless package has been opened or damaged	LOT	Sterile Lot Number
	One use only, do not reuse	SN	Serial Number
Stored in 50% ethanol	Stored in 50% ethanol		Straight Prosthesis
Imported by	Imported by		Curved Prosthesis
EC REP	Authorised Representative in the European Union		Use by
	Internal diameter		Manufactured by
	Minimum length	#	Quantity

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